

# USER MANUAL

## InviTrap<sup>®</sup> Spin Plant RNA Mini Kit

for purification of total RNA from plant material and cells (up to 100 mg or  $1 \times 10^7$  cells)  
and from filamentous fungi

## Instruction for the InviTrap® Spin Plant RNA Mini Kit

The **InviTrap® Spin Plant RNA Mini Kit** is qualified for the isolation of high quality total RNA from plant cells or plant tissues and filamentous fungi (up to 100 mg or  $1 \times 10^7$  cells) in a spin filter format. Special buffer conditions guarantee an efficient lysis of the starting material and a simultaneous inactivation of endogenous RNases. Genomic DNA is separated from the total RNA by binding to specially optimized mineral carrier particles during the lysis. These particles are included in the Lysis Buffer. Multiple samples can be processed simultaneously in less than 1 hour. The kit can be used further for simultaneous isolation of total RNA and proteins.

DNase-digestion to remove contaminating genomic DNA is not necessary.

The **InviTrap® Spin Plant RNA Mini Kit** is intended for life science research use only. Prior to using it for other purposes, the user must validate the system in compliance with the applicable law, directives, and regulations.

The kit is neither validated for the isolation of total RNA from serum, plasma, blood, tissue, human or animal cell culture, bacteria or yeast cells nor for viruses. The performance of the kit in isolating and purifying total RNA from fecal samples has not been evaluated.

**Trademarks:** InviTrap®, Eppendorf®. Registered marks, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

The Invisorb® technology is covered by patents and patent applications: US 6,110,363, US 6,043,354, US 6,037,465, EP 0880535, WO 9728171, WO 9534569, EP 0765335, DE 19506887, DE 10041825.2, WO 0034463.

InviTrap® and Invisorb® are registered trademarks of Invitex Molecular GmbH.

The PCR process is covered by US Patents 4,683,195, and 4,683,202 and foreign equivalents owned by Hoffmann-La Roche AG.

© 2020 Invitex Molecular, all rights reserved.









## Contents

<b>Kit contents of the InviTrap® Spin Plant RNA Mini Kit</b>	<b>4</b>
<b>Symbols</b>	<b>5</b>
<b>Storage</b>	<b>5</b>
<b>Quality control</b>	<b>5</b>
<b>Intended use</b>	<b>6</b>
<b>Product use Limitation</b>	<b>6</b>
<b>Safety information</b>	<b>7</b>
<b>Product characteristics of the InviTrap® Spin Plant RNA Mini Kit</b>	<b>8</b>
<b>Principle and procedure</b>	<b>10</b>
<b>Sampling and storage of starting material</b>	<b>10</b>
<b>Sample homogenization</b>	<b>10</b>
<b>Procedure</b>	<b>11</b>
Lysis	11
Binding and removal of DNA	12
Binding total RNA	12
Removing residual contaminants	12
Elution	12
<b>Important notes</b>	<b>13</b>
Important points before starting a protocol	13
Preparing reagents and buffers	13
Equipment and reagents to be supplied by user	15
<b>Scheme</b>	<b>16</b>
<b>Instructions</b>	<b>17</b>
<i>Protocol 1: Total RNA extraction from plant cells (1x10<sup>7</sup>) or plant tissues and filamentous fungi (max. 100 mg)</i>	17
Supplemental protocol for DNA digestion on the RTA Spin Filter	19
Supplemental Protocol for simultaneous isolation of RNA and proteins	20
<b>Troubleshooting</b>	<b>21</b>
<b>Appendix</b>	<b>23</b>
General notes on handling RNA	23
Storage of RNA	24
Quantification of RNA	24
Purity of RNA	24
<b>Ordering information</b>	<b>27</b>

## Kit contents of the InviTrap® Spin Plant RNA Mini Kit

	50 total RNA extractions	250 total RNA extractions
<b>Catalogue No.</b>	1064100300	1064100400
<b>Lysis Solution DCT</b>	50 ml	250 ml
<b>Lysis Solution RP</b>	50 ml	250 ml
<b>Wash Buffer R1</b>	20 ml (final volume 40 ml)	80 ml (final volume 160 ml)
<b>Wash Buffer R2</b>	2 x 12 ml (final volume 2 x 60 ml)	2 x 50 ml (final volume 2 x 250 ml)
<b>Elution Buffer R</b>	2 x 2 ml	30 ml
<b>Prefilter</b>	50	5 x 50
<b>RTA Spin Filter Set</b>	50	5 x 50
<b>2.0 ml Receiver Tubes</b>	50	5 x 50
<b>Elution Tubes</b>	50	5 x 50
<b>Manual</b>	1	1
<b>Initial steps</b>	<p>Add 20 ml 96-100% ethanol to the bottle <b>Wash Buffer R1</b>.</p> <p>Add 48 ml 96-100% ethanol to each bottle <b>Wash Buffer R2</b>.</p>	<p>Add 80 ml 96-100% ethanol to the bottle <b>Wash Buffer R1</b>.</p> <p>Add 200 ml 96-100% ethanol to each bottle <b>Wash Buffer R2</b>.</p>
	<p><b>Adjust Lysis Solution DCT or Lysis Solution RP with 1/100 volume of 1M DTT or <math>\beta</math>-Mercaptoethanol.</b> Due to the instability of dissolved DTT or <math>\beta</math>-Mercaptoethanol, under oxidative conditions do not mix the whole <b>Lysis Solutions</b> with DTT or <math>\beta</math>-Mercaptoethanol, in case of the kits with 50 and 250 preparations we recommend the preparation of a volume DTT- or <math>\beta</math>-Mercaptoethanol. Containing <b>Lysis Solution</b> shortly before carrying out the purifications adapted to the number of samples that will be processed. Store the remaining <b>Lysis Solution DCT/ Lysis Solution RP</b> and DTT or <math>\beta</math>-Mercaptoethanol separately in accordance to the storage instructions (see "Storage", page 5).</p> <p>It is possible to replace DTT by <math>\beta</math>-Mercaptoethanol. In that case adjust <b>Lysis Solution DCT/ Lysis Solution RP</b> with 1/100 volume of <math>\beta</math>-Mercaptoethanol as described above.</p>	

## Symbols

	Manufacturer
	Lot number
	Catalogue number
	Expiry date
	Consult operating instructions
	Temperature limitation
	Do not reuse
	Humidity limitation

**Attention:** Do not combine components of different kits, unless the lot numbers are identical!

## Storage

All buffers and kit contents of the **InviTrap® Spin Plant RNA Mini Kit**, except **1 M DTT solution** (not provided in the kit) should be stored at room temperature and are stable for at least 12 months.

**Room temperature (RT) is defined as range from 15-30°C.**

Before every use, make sure that all components have room temperature. If there are any precipitates within the provided solutions solve these precipitates by warming carefully (up to 30°C).

**Wash Buffers** charged with ethanol should be appropriately sealed and stored at room temperature.

Store the **1 M DTT solution** at - 20°C to prevent oxidative damage. Under this condition, the solution of **1 M DTT** is stable for 12 months. If the kit components are consumed in more than one run it is recommended to aliquot the **1 M DTT solution** and minimize the number of freezing and thawing cycles.

## Quality control and product warranty

Invitek Molecular warrants the correct function of the **InviTrap® Spin Plant RNA Mini Kit** for applications as described in this manual. Purchaser must determine the suitability of the Product for its particular use. Should any Product fail to perform the applications as described in the manual, Invitek Molecular will check the lot and if Invitek Molecular investigates a problem in the lot, Invitek Molecular will replace the Product free of charge.

Invitek Molecular reserves the right to change, alter, or modify any Product to enhance its performance and design at any time.

In accordance with Invitek Molecular's EN ISO 13485 certified Quality Management System the performance of all components of the **InviTrap® Spin Plant RNA Mini Kit** have been tested separately against predetermined specifications routinely on lot-to-lot to ensure consistent product quality.

If you have any questions or problems regarding any aspects of **InviTrap® Spin Plant RNA Mini Kit** or other Invitek Molecular products, please do not hesitate to contact us. A copy of Invitek Molecular's terms and conditions can be obtained upon request or are presented at the Invitek Molecular webpage.

**For technical support or further information, please contact:**

**from Germany: +49-(0)30-9489-2901/ 2910**

**from abroad: +49-(0)30-9489-2907**

**or contact your local distributor.**

## Intended use

The **InviTrap® Spin Plant RNA Mini Kit** is the ideal tool for reliable and fast manual isolation and purification of high quality total RNA from up to 100 mg or  $1 \times 10^7$  cells of plant material or filamentous fungi. For reproducible high yields an appropriate sample storage and quick operation under the rules for RNA operation is essential. The purified RNA is ready to use.

The isolation protocol and all buffers are optimized to assure a high yield as well as a high purity of purified total RNA. All manual work is reduced to a minimum.

With the **InviTrap®** procedure, all RNA molecules longer than 200 nucleotides are purified quantitatively.

For some studies, a simultaneous investigation of the RNA and of the cell specific proteins is necessary. The kit includes an isolation protocol of RNA and intracellular proteins from the same sample.

**THE PRODUCT IS INTENDED FOR USE BY PROFESSIONALS ONLY, SUCH AS TECHNICIANS, PHYSICIANS AND BIOLOGISTS TRAINED IN MOLECULAR BIOLOGICAL TECHNIQUES.** It is designed to be used with any downstream application employing enzymatic amplification or other enzymatic modifications of RNA followed by signal detection or amplification. Any diagnostic results generated by using the sample preparation procedure in conjunction with any downstream diagnostic assay should be interpreted with regard to other clinical or laboratory findings.

To minimize irregularities in diagnostic results, adequate controls for downstream applications should be used.

## Product use limitation

The kit is neither validated for the isolation of total RNA from serum, plasma, blood, human or animal cell culture or tissue, bacteria or yeast cells, nor for viruses. The performance of the kit in isolating and purifying total RNA from fecal samples has not been evaluated.

The kit was not tested on its ability to desalinate RNA or for RNA purification from enzymatic reactions, like Proteinase digestion, RNA ligation or labeling reactions.

The included chemicals are only useable once.

Differing of starting material or flow trace may lead to inoperability; therefore, neither a warranty nor guarantee in this case will be given, neither implied nor express.

The user is responsible to validate the performance of the Invitek Molecular Product for any particular use. Invitek Molecular does not provide for validation of performance characteristics of the Product with respect to specific applications.

Invitek Molecular products may be used e.g. in clinical diagnostic laboratory systems under following conditions:

- If used in the US, based on the condition that the complete diagnostic system of the laboratory has been validated pursuant to CLIA' 88 regulations.
- For other countries based on the condition that the laboratory has been validated pursuant to equivalents according to the respective legal basis.

All Products sold by Invitek Molecular are subject to extensive quality control procedures (according to EN ISO 13485) and are warranted to perform as described herein. Any problems, incidents or defects shall be reported to Invitek Molecular immediately upon detection thereof.

The chemicals and the plastic parts are for laboratory use only; they must be stored in the laboratory and must not be used for purposes other than intended.

The Product with its contents is unfit for consumption

## Safety information

When and while working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles!

Avoid skin contact! Adhere to the legal requirements for working with biological material!

For more information, please consult the appropriate material safety data sheets (MSDS). These are available online in convenient and compact PDF format at [www.invitek-molecular.com](http://www.invitek-molecular.com) for each Invitek Molecular Product and its components. If buffer bottles are damaged or leaking, **WEAR GLOVES, AND PROTECTIVE GOGGLES** when discarding the bottles in order to avoid any injuries.

Invitek Molecular has not tested the liquid waste generated by the **InviTrap® Spin Plant RNA Mini Kit** procedures for residual infectious materials. Contamination of the liquid waste with residual infectious materials is highly unlikely, but cannot be excluded completely. Therefore, liquid waste must be considered infectious and be handled and discarded according to local safety regulations.

European Community risk and safety phrases for the components of the **InviTrap® Spin Plant RNA Mini Kit** to which they apply, are listed below as follows:

### Lysis Solution DCT



Warning

H302-H319-H332-H412-P280-P305+P351+P338-EUH032

### Lysis Solution RP



Warning

H302-H315-H319-P280-P305+P351+P338

### Wash Buffer R1



Warning

H302-H332-H412-P280-P305+P351+P338-EUH032

H302: Harmful if swallowed.

H315: Causes skin irritation.

H319: Causes serious eye irritation.

H332: Harmful if inhaled.

H412: Harmful to aquatic life with long lasting effects.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

EUH032: Contact with acids liberates very toxic gas.

**Emergency medical information can be obtained 24 hours a day from infotrac:**

**outside of USA: 1 – 352 – 323 – 3500**

**in USA : 1 – 800 – 535 – 5053**

## Product characteristics of the InviTrap® Spin Plant RNA Mini Kit

Starting material	Yield	Time	Ratio
up to 100 mg or 1x10 <sup>7</sup> cells of plant material or filamentous fungi	up to 80 µg total RNA from 100 mg of <i>Oryza sativa</i> cotyledons	about 15 minutes after lysis of the starting material	A <sub>260</sub> : A <sub>280</sub> 1.6 - 2.0

The **InviTrap®** procedure represents a well-established technology for RNA purification. This technology combines the selective binding properties of a silica-based membrane with the speed of microspin technology.

The **InviTrap® Spin Plant RNA Mini Kit** provides a fast and efficient way for reliable isolation of high quality total RNA from up to 100 mg or 1x10<sup>7</sup> cells of plant material or filamentous fungi. Up to 100 mg of sample is homogenized in liquid nitrogen and then lysed under highly denaturing conditions. The provided lysis buffers (**Lysis Solution DCT/ Lysis Solution RP**) of the **InviTrap® Spin Plant RNA Mini Kit** are adapted to the special requirements for the lysis of plant material from different species or different parts of the same plant. Plant material shows a high diversity in its composition and the content of secondary metabolites.

The total RNA isolation combines efficient lysis and homogenization of the starting material in the presence of a highly denaturing guanidiniethiocyanate containing buffer, which immediately inactivates RNases to ensure purification of intact RNA with the efficient removal of the genomic DNA. The genomic DNA is fixed at the surface of the nucleic acid binding mineral carrier particles within the **Lysis Solution DCT** or **Lysis Solution RP** simultaneously to the cell lysis.

Selective binding of the DNA is assured by optimized buffer conditions that reduce the undesired binding of the RNA to a minimum. Some secondary metabolites of plants are able to bind to the DNA-binding carrier particles and block the interaction with the DNA. Thus, the complete DNA removal can be disturbed. To avoid this effect **Lysis Solution DCT** is optimized for plant material showing a high polysaccharide concentration and **Lysis Solution RP** is adapted to plant material with high phenol content (see page 13). Depending on the content and type of secondary metabolites in some plant materials or mycelia of filamentous fungi guanidine thiocyanate can cause solidification of the sample. Therefore, the use of **Lysis Buffer RP** is recommended.

During the lysis procedure under permanent shaking the carrier material is also shearing high molecular components, homogenizing the plant material and reducing the viscosity. After lysis, the mineral carrier particles loaded with the genomic DNA are removed by centrifugation.

Ethanol is added to provide appropriate binding conditions, and the sample is then applied to an RTA Spin Filter where the total RNA binds to the membrane and contaminants are efficiently washed away by repeated wash steps. The purified total RNA can be eluted in a small volume of RNase free water. The isolated total RNA is ready to use and should be stored at - 80°C.

If the subsequent RNA applications show a high sensitivity to small DNA contaminations, an optional DNase digestion is recommended

The extracted RNA contains enriched mRNA. The particle size distribution of purified RNA is similar to those, gained using a CsCl - gradient.



Due to the high purity, the isolated total RNA is ready to use for a broad panel of downstream applications like:

- Northern Blotting,
- RNA dot blots,
- cDNA-library,
- *in vitro* translation,
- RT-PCR\*,
- DDRT-PCR\*,
- or TaqMan® analysis and array technologies.

**For technical support or further information, please contact:**

+49-(0)30-9489-2901/ 2907/ 2910

or your local distributor.

\* The PCR process is covered by US Patents 4,683,195, and 4,683,202 and foreign equivalents owned by Hoffmann-La Roche AG.

## Principle and procedure

The InviTrap® Spin Plant RNA Mini Kit procedure comprises the following steps:

- disruption and homogenization of the plant material
- cell lysis
- selective binding of the genomic DNA to a specific carrier and separation of contaminating DNA on the surface of the DNA binding mineral carrier particles
- transfer of the remaining sample into the RNA binding RTA Spin Filter, followed by the adjustment of the binding conditions
- binding of the total RNA to the membrane, while proteins remain in the filtrate
- washing of the membrane and elimination of contaminants and ethanol
- elution of highly pure total RNA

Repeated wash steps make sure that contaminations and enzyme inhibitors are efficiently removed and high purified RNA is eluted in **Elution Buffer R** or RNase free water.

This manual contains 1 protocol (+ 2 supplemental for protocols) (17-20).

### Sampling and storage of starting material

Best results are obtained using fresh plant material. As long as samples are not shock frosted or are incubated with RNase inhibitors or denaturing reagents, the RNA is not secured. Therefore, it is essential, that cells are immediately flash frozen after cell harvest and are stored at - 80°C. RNA from deep frozen samples is stable for months. Frozen tissue should not be allowed to thaw while handling. RNA purification should be processed as soon as possible. Cells can be stored also in **Lysis Solution DCT/ Lysis Solution RP** at - 80°C after cell lysis. To avoid overloading the kit, weighting tissues or direct counting of cells is recommended to determine the amount of starting material.

Invitex Molecular will be released of its responsibilities if other sample materials than described in the Intended Use are processed or if the sample preparation protocols are changed or modified.

### Sample homogenization

The efficient and complete disruption and homogenization of the plant tissue samples is absolute essential for isolation of a high yield of total RNA.

### Disrupting and homogenizing of the plant material are 2 distinct steps

The **disruption procedure**, the breakage of intercellular matrix, like cells walls, organelles and plasma membranes, is necessary to release the nucleic acids contained in the cell, thus inefficient disruption decreases the RNA yield. Different samples require different methods to achieve complete disruption

The **homogenization** means the reduction of the viscosity of the lysate after disruption. Contaminating genomic DNA and other cellular components of high molecular weight are sheared to form a homogenous lysate. Incomplete homogenization results in inefficient binding of RNA to the RTA Spin Filter membrane and therefore significantly reduced RNA yields.

It is possible to use a commercially available bead mills in combination with or without beads for the disruption and homogenization of the starting material. Alternatively, the starting material can be reduced to a fine powder in liquid nitrogen using a mortar and pestle.

Some disruption methods simultaneously homogenize the sample, while others require an additional homogenization step.

## Disruption and homogenization using the bead-milling system

In bead-milling plant tissues can be disrupted by rapid agitation in the presence of stainless steel or tungsten carbide beads. Disruption and simultaneous homogenization occur by the shearing and crushing action of the beads as they collide with the plant cells. In this case, plant material, beads, and disruption vessels must all be flash frozen in liquid nitrogen, and disruption is performed without lysis buffer

Disruption efficiency is influenced by:

- size and composition of beads
- ratio of buffer to beads
- amount of starting material
- speed and configuration of the bead mill
- disintegration time

## Disruption using a mortar and pestle

For disruption using a mortar and pestle, freeze the plant tissue immediately in liquid nitrogen and grind to a fine powder under liquid nitrogen. Transfer the suspension (plant tissue powder and liquid nitrogen) into a liquid nitrogen cooled, appropriately sized tube and allow the liquid nitrogen to evaporate without allowing the sample to thaw.

Add lysis buffer and continue as quickly as possible with the homogenization according to one of the two methods below.

**Note:** *Grinding the sample using a mortar and pestle will disrupt the sample, but will not homogenize it. Homogenization must be performed afterwards.*

**This homogenization step will be supported in all InviTrap® RNA Spin kits via the mineral carrier particles provided in the Lysis buffers during the lysis step.**

## Determining the correct amount of starting material

It is essential to use the correct amount of starting material in order to obtain optimal RNA yield and purity. A maximum amount of 100 mg plant material or  $1 \times 10^7$  cells can generally be processed. For most plant materials, the RNA binding capacity of the RTA Spin Filter and the lysing capacity of **Lysis Solution DCT/ Lysis Solution RP** will not be exceeded by these amounts.

If there is no information about the nature of your starting material, it is recommended to start with no more than 50 mg plant material or  $3\text{--}4 \times 10^6$  cells. Depending on RNA yield and purity, it may be possible to use up to 100 mg plant material or up to  $1 \times 10^7$  cells in subsequent preparations.

**Do not overload the RTA spin column, as this will significantly reduce RNA yield and quality.**

Counting cells or weighing tissue is the most accurate way to quantitate the amount of starting material. As a guide, a 1.5 cm diameter leaf disc weighs 25–75 mg.

## Procedure

### Lysis

After complete disruption and homogenization of the starting material, the cells are directly charged with the DTT- or  $\beta$ -Mercaptoethanol containing **Lysis Solution DCT** or **Lysis Solution RP** and thoroughly mixed. Shake the **Lysis Solution DCT / Lysis Solution RP** to resuspend the mineral carrier particles prior to use. Due to the strong denaturing lysis conditions cells are quickly broken and RNases are inactivated simultaneously. The RNA is secured. **DTT** or  $\beta$ -**Mercaptoethanol** is added to inactivate the RNases by cleaving intramolecular disulfide bridges.

### **Binding and removal of DNA**

In general, DNase digestion is not necessary to remove the genomic DNA. Under lysis conditions the DNA is efficiently bound by the mineral carrier particles contained in **Lysis Solution DCT/ Lysis Solution RP**. The particle-bound DNA is removed from the solution by subsequent centrifugation. However, RNA applications that show an increased sensitivity to contaminations of DNA (e.g. caused by overloading the kit) require an optional DNase digestion. If the DNase digestion is carried out while running RNA purification protocol the DNase will be removed completely in the wash steps of the RNA binding RTA Spin Filter.

### **Binding total RNA**

After removal of the genomic DNA from the lysate by removal of the DNA binding mineral carrier particles the total RNA is found in the supernatant. To adjust the RNA binding conditions 96 - 100% ethanol is added. Following carefully mixing the solution is completely transferred onto the RTA Spin Filter and the RNA is bound on the membrane by centrifugation.

### **Removing residual contaminants**

Contaminants are efficiently washed away using **Wash Buffer R1** and **R2**, while the RNA remains bound to the membrane.

### **Elution**

Total RNA is eluted from the membrane using 30 - 60 µl **Elution Buffer R** (or RNase free water). The eluted RNA is ready for use in different subsequent applications.

## Important notes

### Important points before starting a protocol

Immediately upon receipt of the Product, inspect the Product and its components as well as the package for any apparent damages, correct quantities and quality. If there are any unconformities you have to notify Invitex Molecular in writing with immediate effect upon inspection thereof. If buffer bottles are damaged, contact the Invitex Molecular Technical Services or your local distributor. In case of liquid spillage, refer to "Safety Information" (see page 7). Do not use damaged kit components, since their use may lead to poor kit performance.

- Always change pipet tips between liquid transfers. To avoid cross-contaminations, we recommend the use of aerosol-barrier pipet tips.
- All centrifugation steps are carried out at room temperature.
- When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles.
- Discard gloves if they become contaminated.
- Do not combine components of different kits unless the lot numbers are identical.
- Avoid microbial contamination of the kit reagents.
- To minimize the risk of infections from potentially infectious material, we recommend working under laminar air-flow until the samples are lysed.
- Trained personnel should only use this kit.

## Preparing reagents and buffers

### 1. Preparing buffers

#### 50 total RNA-extractions:

Adjust **Lysis Solution DCT/ Lysis Solution RP** with 1/100 volume of 1 M DTT\* or  $\beta$ -**Mercaptoethanol**.

add 20 ml 96-100% ethanol to the bottle **Wash Buffer R1**

add 48 ml 96-100% ethanol to each bottle **Wash Buffer R2**

#### 250 total RNA-extractions:

Adjust **Lysis Solution DCT/ Lysis Solution RP** with 1/100 volume of 1 M DTT\* or  $\beta$ -**Mercaptoethanol**.

add 80 ml 96-100% ethanol to the bottle **Wash Buffer R1**

add 200 ml 96-100% ethanol to each bottle **Wash Buffer R2**

\*It is possible to replace DTT by  $\beta$ -Mercaptoethanol. Adjust **Lysis Solution DCT/ Lysis Solution RP** with 1/100 volume of  $\beta$ -Mercaptoethanol.

## 2. DTT (not provided)

Adjust **Lysis Solution DCT** with 1/100 volume of 1 M **DTT** or **β-Mercaptoethanol**. Due to the instability of dissolved **DTT** or **β-Mercaptoethanol** under oxidative conditions do not mix the whole **Lysis Solution R** with **DTT** or **β-Mercaptoethanol** in case of the kits with 50 and 250 preparations. We recommend the preparation of a volume **DTT**-containing **Lysis Solution R** shortly before carrying out the purifications adapted to the number of samples that will be processed. Store the remaining **Lysis Solution DCT** and **DTT** or **β-Mercaptoethanol** separately in accordance to the storage instructions (see “Storage”, page 5).

It is possible to replace **DTT** by **β-Mercaptoethanol**. In that case, adjust **Lysis Solution DCT** with 1/100 volume of **β-Mercaptoethanol** as described above.

## 3. Lysis Solution

Shake **Lysis Solution DCT/ Lysis Solution RP** gently before use to homogenize the DNA-binding mineral carrier particles! Wait a short time because of foam formation! Cell lysates stored in **Lysis Solution DCT/ Lysis Solution RP** are stable at - 80°C for several months.

### Recommendations for the use of

<b>Lysis Solution RP</b>	<b>Lysis Solution DCT</b>
for phenol containing plants	for polysaccharide containing plants
<i>Allium cepa</i> <sup>3</sup>	<i>Actinidia deliciosa</i>
<i>Arabidopsis thaliana</i> <sup>1</sup>	<i>Beta vulgaris</i>
Ascomyceten	<i>Brassica napus</i> <sup>1</sup>
<i>Brassica oleracea</i> <sup>1</sup>	Chlorophyceae
<i>Casuarina</i> <sup>1/3/5</sup> (Australien tree)	Fungi
<i>Fusarium avenaceum</i> <sup>4</sup>	<i>Hordeum vulgare</i> <sup>1</sup>
<i>Gossypium spec.</i> <sup>1</sup>	<i>Solanum tuberosum</i> <sup>1</sup>
<i>Lycopersicon esculentum</i> <sup>1</sup>	<i>Arabidopsis thaliana</i> <sup>1</sup>
<i>Malus sylvestris</i> <sup>1</sup>	<i>Nicotiana tabacum</i> BY-2
<i>Mangifera indica</i>	
<i>Nicotiana tabacum</i> <sup>1</sup>	
<i>Oryza sativa</i> <sup>2</sup>	
<i>Picea abies</i>	
<i>Rosa hybrida</i> <sup>1</sup>	
Solanaceae <sup>1</sup>	
<i>Vitis vinifera</i> <sup>1</sup>	
<i>Zea mays</i>	

<sup>1</sup> Leave; <sup>2</sup> Cotyledone; <sup>3</sup> Root; <sup>4</sup> Mycelium; <sup>5</sup> Bark

## Equipment and reagents to be supplied by user

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information, please consult the appropriate material safety data sheets (MSDS). (See our webpage: [www.invitek-molecular.com](http://www.invitek-molecular.com))

- Microcentrifuge ( $\geq 12.000$  rpm)
- Ethanol (96-100%)
- Sterile; RNase-free pipet tips
- Disposable gloves
- DTT or  $\beta$ -Mercaptoethanol.

Possible suppliers:	Eppendorf AG 22331 Hamburg, Germany Phone: +49 (0) 40 53801 0 Fax: +49 (0) 40 53801 556 E-Mail: <a href="mailto:eppendorf@eppendorf.com">eppendorf@eppendorf.com</a> Internet: <a href="http://www.eppendorf.com">www.eppendorf.com</a>	SIGMA Laborzentrifugen GmbH 37507 Osterode am Harz, Germany Phone: +49-5522-5007-0 Fax: +49-5522-5007-12 E-Mail: <a href="mailto:info@sigma-zentrifugen.de">info@sigma-zentrifugen.de</a> Internet: <a href="http://www.sigma-zentrifugen.de">www.sigma-zentrifugen.de</a>
---------------------	--	---

### **Bead-mill homogenizers & stainless steel and tungsten carbide beads can be purchased from:**

QIAGEN ([www.qiagen.com](http://www.qiagen.com)) TissueLyser system  
UniEquip ([www.uniequip.com](http://www.uniequip.com))

### **Glass, stainless steel, and tungsten carbide beads can be purchased from:**

Retsch ([www.retsch.de](http://www.retsch.de))

### **Rotor-stator homogenizers can be purchased from:**

BioSpec Products, Inc. ([www.biospec.com](http://www.biospec.com)): Tissue-Tearor™ homogenizer  
Charles Ross & Son Company ([www.mixers.com](http://www.mixers.com))  
IKA ([www.ika.de](http://www.ika.de)): ULTRA-TURRAX® dispersers  
KINEMATICA AG ([www.kinematica.ch](http://www.kinematica.ch)) or Brinkmann Instruments, Inc. ([www.brinkmann.com](http://www.brinkmann.com)): POLYTRON® laboratory dispersing devices  
Omni International, Inc. ([www.omni-inc.com](http://www.omni-inc.com))  
Silverson ([www.silverson.com](http://www.silverson.com))  
VirTis ([www.virtis.com](http://www.virtis.com))





## Instructions

The following notes are valid for all protocols:

**Note:** *The centrifugation steps were made with the **Centrifuge 5415 D** from Eppendorf. The indicated **rpm amounts** are referring to this centrifuge.*

### **Protocol 1: Total RNA extraction from plant cells ( $1 \times 10^7$ ) or plant tissues and filamentous fungi (max. 100 mg)**

---

#### **1. Cell/ tissue disruption**

The complete disruption and homogenization of the plant tissue samples is absolute essential for isolation of a high yield of total RNA. The disruption procedure, the breakage of cell walls and plasma membranes, is necessary to release the nucleic acids contained in the cell, thus inefficient disruption decreases the RNA yield. The homogenization means the reduction of the viscosity of the lysate after disruption. Contaminating genomic DNA and other cellular components of high molecular weight are sheared to form a homogenous lysate. If the homogenization of the starting material is not done very carefully the yield of total RNA purified is reduced significantly.

Grind sample under liquid nitrogen thoroughly to a fine powder using a mortar and pestle or use bead mills. Incomplete grinding of the plant material will lead to reduced total RNA yields! Transfer the tissue powder together with liquid nitrogen to an appropriately sized tube (e.g. 2 ml reaction tube) and allow the liquid nitrogen to evaporate. Continue immediately with step 2. To avoid total RNA-degradation, do not allow the sample to thaw.

#### **2. Cell Lysis**

***Important:*** *Shake Lysis Solution DCT or Lysis Solution RP gently before use! Add DTT or  $\beta$ -Mercaptoethanol!!! Wait a short time because of foam formation!*

Add 900  $\mu$ l DTT or  $\beta$ -Mercaptoethanol containing **Lysis Solution DCT** or **Lysis Solution RP** to a maximum of 100 mg tissue powder. Mix thoroughly by vortexing and incubate for 15 - 30 min at room temperature in a thermomixer under continuously shaking or vortexing from time to time. (Best results will be achieved, if the sample are incubated under shaking conditions !) In some applications a short (2-3 min) incubation step at 55°C may lead to higher total RNA yields. However for samples with high starch concentrations, incubation at high temperatures should be omitted to prevent swelling of the starting material.

#### **3. Removing of the genomic DNA**

Centrifuge the sample at maximum speed for 1 min. After centrifugation a pellet is formed. Place the Prefilter into a 2.0 ml Receiver Tube. Transfer the cleared supernatant carefully into the Prefilter. Centrifuge for 1 min at 11.000 x g (11.000 rpm). Discard the Prefilter. The filtrate (app. 800  $\mu$ l) contains the RNA.

#### **4. Adjust RNA binding conditions**

Add 500  $\mu$ l of 96-100 % **ethanol** to the filtrate. Mix very well by pipetting some times up and down. It is important to mix the lysate completely with the ethanol.

***Important:*** *Other ethanol concentration than 96 –100 % will lead to a significantly reduced yield!*

## 5. Binding of the total RNA to the RTA Spin Filter

Transfer app. 750 µl of the lysate to a RTA Spin Filter. Incubate for 1 min and centrifuge at 11.000 x g (11.000 rpm) for 2 min. Discard the flow-through and place the RTA Spin Filter back into the Receiver Tube. Reload the RTA Spin Filter with the residual volume of lysate and centrifuge again for 1 min. Discard the flow-through and place the RTA Spin Filter back into the RTA Receiver Tube.

## 6. First washing of the RTA Spin Filter

Add 500 µl **Wash Buffer R1** onto the RTA Spin Filter and centrifuge for 1 min at 11.000 x g (10.000 rpm). Discard the flow-through and reuse the RTA Receiver Tube.

## 7. Second washing of the RTA Spin Filter

Add 700 µl **Wash Buffer R2** onto the RTA Spin Filter and centrifuge for 1 min at 11.000 x g (10.000 rpm). Discard the flow-through and reuse the RTA Receiver Tube.

**Repeat this washing step once!**

## 8. Drying of the RTA Spin Filter

To eliminate any traces of ethanol, centrifuge for 4 min at 11.000 x g (11.000 rpm). Discard the RTA Receiver Tube.

## 9. Elution of total RNA

Transfer the RNA RTA Spin Filter into a RNase-free Elution Tube and add 30 - 60 µl of **Elution Buffer R** directly onto the membrane of the RTA Spin Filter. Incubate for 2 min and centrifuge for 1 min at 11.000 x g (11.000 rpm). Discard the RTA Spin Filter and place the eluted total RNA immediately on ice!

**Note:** *Depending on the expected yield or the needed concentration of total RNA you can elute with different volumes of **Elution Buffer R**. A lower volume of **Elution Buffer R** will increase the concentration of RNA and a higher volume of **Elution Buffer R** will lead to an increased yield but a lower concentration of total RNA. Please note, that the minimum of **Elution Buffer R** should be 30 µl.*

## Supplemental protocol for DNA digestion on the RTA Spin Filter

---

If you want to get a totally DNA free sample for RT-qPCR Applications or other protocols, where smallest amounts of DNA might disturb you may follow this protocol and do a DNase digestion directly on the RTA Spin Filter.

**Therefore, you have to change Step 7 in protocol 1 according to the following procedure:**

### **7. Second washing of the RTA Spin Filter**

Add 700  $\mu$ l **Wash Buffer R2** onto the RTA Spin Filter with the bound RNA and centrifuge for 1 min at 11.000 x g (11.000 rpm) Discard the flow-through and put the RTA Spin Filter back into the RTA Receiver Tube.

Add 10  $\mu$ l of a DNase reaction mixture directly in the center of the RTA Spin Filter membrane. A typical reaction mix is for example 1  $\mu$ l of **DNase I** (50 u) enzyme in 9  $\mu$ l 1 x DNase reaction buffer, but here follow the advice of the manufacturer. Do not exceed the reaction mix volume to more than 10  $\mu$ l! Incubate the RTA Spin Filter at room temperature for 10 min. Add again 600  $\mu$ l **Wash Buffer R2** onto the RTA Spin Filter incubate for 1 minute and centrifuge for 1 min sec at 11.000 x g (11.000 rpm). Reuse the Receiver Tube.

**Repeat this washing step once again and follow the protocol further from step 8.**

## Supplemental protocol for simultaneous isolation of RNA and proteins

---

The starting material for RNA purifications is often limited and in some cases the isolation of proteins from the same sample is necessary. Following the given protocol below the simultaneous preparation of RNA and proteins is possible using the **InviTrap® Spin Plant RNA Mini Kit**.

**Important Note:** *The buffer solutions for RNA purification provide very denaturing conditions due to the contained salts. Proteins isolated according to the protocol below are denatured and can be analyzed by SDS-PAGE or Western Blot. Experiments requiring the native state of the proteins (e.g. interaction assays) are excluded.*

### 1. Collection of protein containing material

The protein fraction is found in the flow-through of the RTA Spin Filter in the purification protocol (see *Protocol 1*, step 5). Measure the volume of the flow-through.

### 2. Precipitation of proteins

Add the three fold volume of ice cold **acetone** and mix well by vortexing.  
Centrifuge for 10 min at 13.400 x g (12.000 rpm) at 4°C.  
Discard the supernatant. Be careful not to remove the pellet!

### 3. Washing step

Add 500 µl of ice cold **ethanol** and centrifuge for 3 min at 13.400 x g (12.000 rpm) at 4°C.  
Discard the supernatant, be careful not to remove the pellet.

### 4. Dissolving of proteins

Dissolve the protein pellet by suspending in a buffer solution suitable for the subsequent application. For SDS-PAGE, directly dissolve proteins in 1-fold Laemmli Buffer and heat them at 99°C for 5 min.

**ATTENTION:** *Never try a TCA-Precipitation from RNA-Lysis Buffers, poisonous gas will be generated!*

## Troubleshooting

Problem/ Probable cause	Comments and suggestions
<p><b>clogged RTA Spin Filter</b></p> <p>insufficient disruption or homogenization of the starting material</p> <p>too much starting material</p>	<p>after lysis spin lysate to pellet debris and continue with the protocol using the supernatant</p> <p>increase g-force and/ or centrifugation time</p> <p>all centrifugation steps should be conducted at room temperature</p> <p>in subsequent preparations, reduce the amount of starting material and/ or increase volume of <b>Lysis Solution DCT</b> or <b>RP</b></p>
<p>very viscous lysate, like a gel after addition of <b>Lysis Solution DCT</b> or <b>RP</b> to plant material</p>	<p>if too many plant material including DNA has been used, the carrier with the bound DNA will be too viscous to pipet after homogenization</p> <p>in this case, divide the sample into two aliquots and adjust the volumes of each aliquot to 900 µl with <b>Lysis Solution DCT</b> or <b>RP</b>; continue with the procedure from step</p>
<p><b>low A<sub>260</sub>/A<sub>280</sub> value</b></p> <p>RNA sample is diluted in H<sub>2</sub>O</p>	<p>do not use RNase free water to dilute the sample for measuring the RNA purity</p> <p>the use of a neutral buffer (10 mM Tris/HCl, pH 7.0) is recommended</p>
<p>protein contamination of RNA</p>	<p>use less sample</p> <p>repeat washing step one</p>
<p>RNA sample is too strong diluted</p>	<p>sample concentration is out of the range of the photometer</p>
<p><b>little or no total RNA eluted</b></p> <p>incomplete removal of culture medium from fungi culture etc.</p>	<p>make sure that the cell culture medium is complete removed before starting the isolation</p>
<p>insufficient disruption or homogenization</p>	<p>reduce the amount of starting material, overloading the kit reduces the yield</p>
<p>incomplete elution</p>	<p>prolong the incubation time with <b>Elution Buffer R</b> to 5-10 min or repeat elution step once again</p> <p>elute a second time with 100 µl <b>Elution Buffer R</b></p> <p>incubate RTA Spin Filter at room temperature (15-25°C) for 5 min with RNase free water or <b>Elution Buffer R</b> prior centrifugation</p>
<p>no alcohol added to the <b>Wash Buffer R1</b> and <b>R2</b></p>	<p>check that <b>Wash Buffer R1</b> and <b>Wash Buffer R2</b> concentrates were diluted with correct volume of 96-100% ethanol. Do not use denatured alcohol, which contains other substances such as methanol or methylethylketone,</p> <p>repeat the purification procedure with a new sample</p>

<p><b>DNA-contamination</b></p> <p>too much starting material</p> <p>no optimal homogenization of DNA binding carrier</p> <p>plant material contains too many component disturbing DNA binding</p>	<p>reduce amount of starting material. DNase digestion of the eluate containing the total RNA</p> <p>shake <b>Lysis Solution DCT/ Lysis Solution RP</b> carefully before use</p> <p>in subsequent preparations, reduce the amount of starting material and/ or increase volume of <b>Lysis Solution DCT or RP</b></p> <p>actual perform a DNase digestion from the eluate</p>
<p><b>RNA degraded</b></p> <p>inappropriate handling of the starting material</p>	<p>the RNA purification protocol should be performed quickly (see also "General notes on handling RNA", page 23)</p> <p>plant material stored at - 80°C for later processing should be immediately frozen after cell harvest by liquid nitrogen treatment</p>
<p><b>Lysis Solution DCT/ RP does not contain DTT or β-Mercaptoethanol.</b></p>	<p>ensure that <b>DTT</b> or <b>β-Mercaptoethanol</b>.has been added to the <b>Lysis Solution DCT/ Lysis Solution RP</b></p>
<p>RNase contamination</p>	<p>check for RNase contamination of buffers</p> <p>be certain not to introduce any RNase throughout the procedure or during further handling for analysis, see appendix , "General Notes on handling RNA" (page 23)</p>
<p><b>total RNA does not perform well in downstream-applications (e.g. RT-PCR)</b></p> <p>ethanol carryover during elution.</p> <p>salt carryover during elution.</p>	<p>increase g-force or centrifugation time when drying the RTA Spin Filter</p> <p>ensure that <b>Wash Buffer R1</b> and <b>R2</b> are at room temperature</p> <p>check up <b>Wash Buffer R1</b> and <b>R2</b> for salt precipitates. If there are any precipitates solve these precipitates by careful warming</p>

## Appendix

### General notes on handling RNA

RNA is far less stable than DNA. It is very sensitive to degradation by endogenous RNases in the biological material and exogenous RNases, which are permanently present everywhere in the lab. To achieve satisfactory qualitative and quantitative results in RNA preparations, contaminations with exogenous RNases has to be reduced as much as possible. Avoid handling bacterial cultures, cell cultures or other biological sources of RNases in the same lab where the RNA purification is to be carried out.

All glassware should be treated before use to ensure that it is RNase free. Glassware should be cleaned with detergent, thoroughly rinsed and oven baked at 240°C for four or more hours before use. Autoclaving alone will not completely inactivate many RNases. Oven baking will both inactivate RNases and ensure that no other nucleic acids (such as Plasmid DNA) are present on the surface of the glassware. You can also clean glassware with 0.1% DEPC (diethyl pyrocarbonate). The glassware has to stand 12 hours at 37°C and then autoclave or heat to 100°C for 15 min to remove residual DEPC.

- Electrophoresis tanks should be cleaned with detergent solution (e.g. 0.5% SDS), thoroughly rinsed with RNase-free water, and then rinsed with ethanol and allowed to dry.
- Non-disposable plasticware should be treated before use to ensure that it is RNase free. Plastic ware should be thoroughly rinsed with 0.1 M NaOH, 1 mM EDTA followed by RNase free water. You can also take chloroform-resistant plastic ware rinsed with chloroform to inactivate RNases.
- All buffers must be prepared with RNase-free ddH<sub>2</sub>O.
- When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles.
- Change gloves frequently and keep tubes closed.
- All centrifugation steps are carried out at room temperature
- To avoid cross contamination cavity seams should not be moisten with fluid.
- Reduce the preparation time as much as possible.
- Use only sterile, disposable polypropylene tubes throughout the procedure (these tubes are generally RNase free).
- Keep isolated RNA on ice.
- Do not use kit components from other kits with the kit you are currently using, unless the lot numbers are identical.
- To minimize the risk of infections from potentially infectious material, we recommend working under laminar air-flow until the samples are lysed.
- Trained personnel should only use this kit.

## Storage of RNA

Purified viral RNA can be stored at -80°C and is stable for several years at this condition.

## Quantification of RNA

The concentration of RNA should be determined by measuring the absorbance at 260 nm ( $A_{260}$ ) in a spectrophotometer. Readings should be greater than 0.10 to ensure significance. An absorbance of 1 unit at 260 nm corresponds to 40 µg of RNA per ml. This relation is valid only for measurements at neutral pH. The ratio between absorbance values at 260 nm and 280 nm gives an estimate of RNA purity (see below).

When measuring RNA samples, make sure that cuvettes are RNase-free, esp. if the RNA is to be recovered after spectrophotometry. This can be accomplished by washing cuvettes with 0.1 N NaOH, 1 mM EDTA followed by washing with RNase-free water. Use the buffer in which the RNA is diluted for calibration of the spectrophotometer.

An example of the calculation involved in RNA quantification:

- volume of RNA sample: 100 µl
- dilution = 20 µl of RNA sample + 180 µl of 10 mM Tris/HCl pH 7.0 (1/10 dilution).
- measure absorbance of diluted sample in a 0.2 ml cuvette (RNase free):  $A_{260} = 0.2$

Concentration of the RNA sample           = 40 µg/ml \*  $A_{260}$  \* dilution factor  
  = 40 µg/ml \* 0.2 \* 10  
  = 80 µg/ml

Total amount                                   = concentration \* volume of sample in ml  
  = 80 µg/ml \* 0.1 ml  
  = 8 µg of RNA

## Purity of RNA

The ratio of the readings at 260 nm and 280 nm ( $A_{260}/A_{280}$ ) provides an estimate of the purity of RNA with respect to the contaminants that absorb in the UV, such as protein. However, the  $A_{260}/A_{280}$  ratio is influenced considerably by pH. Since water is not buffered, the pH and the resulting  $A_{260}/A_{280}$  ratio can vary greatly. Lower pH results in lower  $A_{260}/A_{280}$  ratio and reduced sensitivity to protein contaminations.\* For accurate values, it is recommended to measure absorbance in 10 mM Tris Cl, pH 7.5. Pure RNA has an  $A_{260}/A_{280}$  ratio of 1.9-2.1\*\* in 10 mM Tris/HCl pH 7.5. Always be sure to calibrate the spectrophotometer with the same solution.

For determination of the RNA concentration, however, it is recommended diluting the sample in a buffer with neutral pH since the relationship between absorbance and concentration ( $A_{260}$  reading of 1 = 40 µg/ml of RNA) is based on an extinction coefficient calculated for RNA at neutral pH.

\* Wilfinger, W.W., Mackey, M., and Chomczynski, P. (1997) Effect of pH and ionic strength on the spectrophotometric assessment of nucleic acid purity. *Bio Techniques* 22, 474

\*\*Values up to 2.3 are routinely obtained for pure RNA (in 10 mM Tris/HCl, pH 7.5) with some photospectrometers.



## Denaturing agarose gel electrophoresis of RNA

The overall quality of an RNA preparation may be assessed by electrophoresis on a denaturing agarose gel; this will also give some information about RNA yield. A denaturing gel system is suggested because most RNA forms extensive secondary structure via intramolecular base pairing, and this prevents it from migrating strictly according to its size. Be sure to include a positive control RNA on the gel so that unusual results can be attributed to a problem with the gel or a problem with the RNA under analysis. RNA molecular weight markers, RNA sample known to be intact, or both, can be used for this purpose.

The given denaturing agarose gel method for RNA electrophoresis is modified from "Current Protocols in Molecular Biology", Section 4.9 (Ausubel et al., eds.).

### 1. Prepare the gel.

- a. heat 1 g agarose in 72 ml water until dissolved, then cool to 60°C.
- b. add 10 ml 10X MOPS running buffer, and 18 ml 37% formaldehyde (12.3 M).

*WARNING: Formaldehyde is toxic through skin contact and inhalation of vapors. Manipulations involving formaldehyde should be done in a chemical fume hood.*

10x MOPS running buffer:   0.4 M MOPS, pH 7.0  
                                  0.1 M sodium acetate  
                                  0.01 M EDTA

- c. pour the gel using a comb that will form wells large enough to accommodate at least 25 µl.
- d. assemble the gel in the tank, and add enough 1X MOPS running buffer to cover the gel by a few millimeters. Then remove the comb.

### 2. Prepare the RNA sample.

**Heat denature samples at 65-70°C for 5-15 min.**

Denaturation for 5 min is typically sufficient for simply assessing RNA on a gel, but a 15 min denaturation is recommended when running RNA for a Northern Blot. The longer incubation may be necessary to completely denature the RNA.

To 1-3 µg RNA, add 0.5-3X volumes Formaldehyde Load Dye.

To simply check the RNA on a denaturing gel, as little as 0.5 X Formaldehyde Load Dye can be used, but to completely denature the RNA, e.g. for Northern blots, use 3 X volumes of Formaldehyde Load Dye. Ethidium bromide can be added to the Formaldehyde Load Dye at a final concentration of 10 µg/ml. Some size markers may require significantly more than 10 µg/ml ethidium bromide for visualization. To accurately size your RNA, however, it is important to use the same amount of ethidium bromide in all the samples (including the size marker) because ethidium bromide concentration affects RNA migration in agarose gels.


### 3. Electrophoresis

Load the gel and electrophorese at 5-6 V/cm until the bromophenol blue (the faster-migrating dye) has migrated at least 2-3 cm into the gel, or as far as 2/3 the length of the gel.

#### 4. Results

Visualize the gel on a UV transilluminator. (If ethidium bromide was not added to the Formaldehyde Load Dye, the gel will have to be post-stained and destained.)

Intact total RNA run on a denaturing gel will have sharp 28S and 18S rRNA bands (eukaryotic samples). The 28S rRNA band should be approximately twice as intense as the 18S rRNA band (Figure 1, lane 3). This 2:1 ratio (28S:18S) is a good indication that the RNA is intact. Partially degraded RNA will have a smeared appearance, will lack the sharp rRNA bands, or will not exhibit a 2:1 ratio. Completely degraded RNA will appear as a very low molecular weight smear (Figure 1, lane 2). Inclusion of RNA size markers on the gel will allow the size of any bands or smears to be determined and will also serve as a good control to ensure the gel was run properly (Figure 1, lane 1). Note: Poly(A) selected samples will not contain strong rRNA bands and will appear as a smear from approximately 6 kb to 0.5 kb (resulting from the population of mRNAs, and depending on exposure times and conditions), with the area between 1.5 and 2 kb being the most intense (this smear is sometimes apparent in total RNA samples as well).



**Figure 1. Intact vs. Degraded RNA.** Two  $\mu\text{g}$  of degraded total RNA and intact total RNA were run beside Ambion's RNA Millennium Markers™ on a 1.5% denaturing agarose gel. The 18S and 28S ribosomal RNA bands are clearly visible in the intact RNA sample. The degraded RNA appears as a lower molecular weight smear.

Generally, at least 200 ng of RNA must be loaded onto a denaturing agarose gel in order to be visualized with ethidium bromide. Some RNA preparations, such as those from needle biopsies or from laser capture microdissected samples, result in very low yields. In these cases, it may be impossible to spare 200 ng of RNA to assess integrity. Alternative nucleic acid stains, such as SYBR® Gold and SYBR® Green II RNA gel stain from Molecular Probes, offer a significant increase in sensitivity over ethidium bromide. Using a 300 nm transilluminator (6 x 15-watt bulbs) and a special filter, as little as 1 ng and 2 ng of RNA can be detected with SYBR® Gold and SYBR® Green II RNA gel stain, respectively.

## Ordering information

<b>Product</b>	<b>Package size</b>	<b>Catalogue No.</b>
InviTrap® Spin Plant RNA Mini Kit	50 preparations	1064100300
InviTrap® Spin Plant RNA Mini Kit	250 preparations	1064100400

## Related products

Invisorb® Spin Plant Mini Kit	250 preparations	1037100300
-------------------------------	------------------	------------

**INVITEK**  
Molecular

Invitek Molecular GmbH  
Röbert-Rössle-Str. 10  
13125 Berlin

Phone: +49 30 94 89 29 01  
Fax: +49 30 94 89 29 09  
[info@invitek-molecular.com](mailto:info@invitek-molecular.com)

[www.invitek-molecular.com](http://www.invitek-molecular.com)

0061064100 V-01-2020